

ePipette

**Intuitive control,
Hold easy,
Work easy!!**



- ✔ Ergonomic and lightweight design. Let you work longer without fatigue.
- ✔ Large, clear OLED-display. Easy setting and easy control.
- ✔ Intuitive user interface that you'll learn how to use it in a few minutes.
- ✔ Always ready. Pipette can be used while recharging with AC adaptor.
- ✔ Fit different brand tip in the market.

AVE  ENE



- Direct charging option with AC-adaptor.
- Large clear OLED-display with all the information you need.
- Intuitive, easy-to-learn interface.
- Reduced force tip ejector.
- Dramatic ergonomic design with minimal effort.
- The autoclavable lower part can easily be removed for decontamination purposes.

ePipette Single Channel Performance Specification

Product Name	Channel	Volume Range (μL)	Increment (μL)	Test Volume (μL)	Inaccuracy %	Imprecision, CV %
ePipette S10	1	0.5 - 10	0.01	1	± 2.5	≅ 1.80
				5	± 1.5	≅ 0.80
				10	± 1.0	≅ 0.40
ePipette S20	1	2 - 20	0.1	2	± 7.5	≅ 2.00
				10	± 1.5	≅ 0.50
				20	± 1.0	≅ 0.30
ePipette S200	1	20 - 200	1	20	± 2.5	≅ 1.00
				100	± 0.8	≅ 0.25
				200	± 0.8	≅ 0.15
ePipette S1000	1	200 - 1000	1	200	± 3.0	≅ 0.60
				500	± 0.8	≅ 0.20
				1000	± 0.8	≅ 0.12

ePipette Ordering Information

Cat No.	Product
AG-PIS0010	ePipette S10
AG-PIS0020	ePipette S20
AG-PIS0200	ePipette S200
AG-PIS1000	ePipette S1000
AG-PIB1100	Spare battery
AG-PICS001	Charging stand (single)

Distributor Info