


# LIAS *Slite*140

*Simple and Smart*



- 1.4 MEGAPIXEL
- MOTORIZED OPTICS
- 8" INCH TOUCH SCREEN
- FRONT & BACK USB PORTS



## SPECIFICATION

CAMERA	
Resolution	1.4Megapixels
Lens	F1.2
Zoom	8.5 – 51, motorized
Exposure control	~10 seconds
Max gel size	16 x 21 cm
Data types	TIFF, JPEG, BMP
Image storage	USB, network

ILLUMINATION	
Trans UV	302 nm
Trans white light	Fold down white light table
Epi white light	Yes

APPLICATIONS	
Fluorescence	Ethidium bromide, SYBR Safe, SYPRO Ruby, etc.
Absorbance/ colorimetry	Silver stain, Coomassie blue, conoly plate

PRODUCT DIMENSIONS	
Footprint	31W x 44L x 66 H cm
Weights	26 kg